

# De Mayo Hip Positioner® Systems

## Accommodates Small to Obese Patients

Building the positioner around a patient means you never need to lift a patient for lateral positioning again. The Bariatric Extension Plate is ideal for obese patients.


Smaller patients are also held securely by the De Mayo Hip Positioner®.

*Secure obese patients without lifting*



US Patent No. 6,820,621

*Foam pads protect patients securely*



(Accommodates the obese patient)

*Security without abdominal pressure*

*Unobstructed access to surgical site*



A single, secure positioning point results in no abdominal pressure for an obese patient.

Crossover bar is tilted away from the hip area, providing clear, unobstructed access to the surgical site and excellent exposure to the acetabulum and femur.

# De Mayo Hip Positioner® Systems

Choose the Complete System that's Best for Your Patient

*"The arm never gets in the way of full access to the acetabulum and femoral canal and also allows you to test range of motion further than any other device we've used."*

*The access is great."*

- Dr. David F. Bindelglass, M.D.

Simple, adjustable positioning for any sized patient on any OR table


The family of complete De Mayo Hip Positioning Systems gives you everything you need to position patients faster, more securely and with less strain on personnel.

Developed by an Orthopedic surgeon, the system is based on a secure patient positioning device that enables you to "build" the positioner around the patient. The collection of stabilizing bars and support pads give you the right option for any patient body type - even your most challenging obese patients. You get unobstructed access to the surgical site with support and security for obese patients that you can't find in any other device.

Latex free single-use padding protects patients and eliminates cross-contamination possibilities. You get a complete system that saves you time and lowers risk for every hip procedure.

## De Mayo Hip Positioner® - LGE


Shown With Bariatric Extension



Support obese patients with no lifting of patient.

## De Mayo Hip Positioner® - STD


Stable Patient Support At The Iliac Crest



Dual compression foam pad provides ultimate patient protection at the iliac crest.



Optional stabilizing bar for extra stability.



imp® Rail Clamp

Large knob provides greater locking ability and attaches to all OR tables.



Complete De Mayo Systems include:

Positioner and all the accessories you need for easier, more secure patient positioning.

Use only imp® Latex Free Single-Use Patient Protective Pads to minimize liability risks while improving positioning stability.

Latex Free



\*Certificate of Liability on file

## De Mayo Hip Positioner® System - LGE

Shown With Optional Stabilizing Bar Secures Obese Anatomies

## De Mayo Hip Positioner® System - STD


Shown With Optional Adjustable Lateral Support Arm With Patient Pad

*"From establishing the position and securing the patient to performing procedures and testing range of motion, we created a system that's faster and easier while also giving doctors completely unrestricted access to the surgical field."*

- Ed De Mayo M.D., device inventor


# De Mayo Hip Positioner® Systems

## Ordering Information

 <b>1101- COMP-SYS</b> <b>NEW COMPLETE SYSTEM</b> <b>DE MAYO HIP POSITIONER®</b>		
QUANTITY	CATALOG #	DESCRIPTION
1	1101-Case	Storage Case
1	1101-BE	Bariatric Extension
1	7516	Locking Bar with Knob
1	1101-A	Lumbar Support
1	1101-B	Crossover Bar (Std)
1	1101-C	Anterior Support
1	1101-D	Anterior Vertical Support (Std)
1	1101-E	Anterior Pelvic Support
1	1101-L-D	Anterior Vertical Support (Lge)
1	1101-L-B	Crossover Bar (Lge)
1	1101-SB-1	Stabilizing Bar (Tube)
1	1101-SB-2	Stabilizing Bar (Hook)
1	707-L	Adjustable Lateral Support w/11" Arm
3	409	<i>imp</i> ® Rail Clamp

<b>1101-L SYS KIT</b> <b>OBESE PATIENT SYSTEM/STABILIZING BAR</b>		
QUANTITY	CATALOG #	DESCRIPTION
1	1101-Case	Storage Case
1	1101-BE	Bariatric Extension
1	7516	Locking Bar with Knob
1	1101-A	Lumbar Support
1	1101-C	Anterior Support
1	1101-E	Anterior Pelvic Support
1	1101-L-D	Anterior Vertical Support (Lge)
1	1101-L-B	Crossover Bar (Lge)
1	1101-SB-1	Stabilizing Bar (Tube)
1	1101-SB-2	Stabilizing Bar (Hook)
2	409	<i>imp</i> ® Rail Clamp

<b>PATIENT PROTECTIVE PADS</b>		
QUANTITY	CATALOG #	DESCRIPTION
10/Case	1101-P	De Mayo Hip Positioner
12/Case	707-P	Adjustable Lateral Support

 <b>1101- SYS KIT</b> <b>STANDARD W/ADJ LATERAL SUPPORT</b>		
QUANTITY	CATALOG #	DESCRIPTION
1	1101-Case	Storage Case
1	7516	Locking Bar with Knob
1	1101-A	Lumbar Support
1	1101-B	Crossover Bar (Std)
1	1101-C	Anterior Support
1	1101-D	Anterior Vertical Support (Std)
1	1101-E	Anterior Pelvic Support
1	707-L	Adjustable Lateral Support w/11" Arm
2	409	<i>imp</i> ® Rail Clamp

<b>1101 KIT</b> <b>DE MAYO HIP POSITIONER - STD ONLY</b>		
---	--	--

<b>1101-L KIT</b> <b>DE MAYO HIP POSITIONER - LGE ONLY</b>		
---	--	--

<b>1101-UPGRADE SYSTEM KIT</b> <b>UPGRADE SYSTEM-CUST WITH 1101-1 OR 1101</b>		
QUANTITY	CATALOG #	DESCRIPTION
1	1101-Case	Storage Case
1	1101-BE	Bariatric Extension
1	1101-L-D	Anterior Vertical Support (Lge)
1	1101-L-B	Crossover Bar (Lge)
1	1101-SB-1	Stabilizing Bar (Tube)
1	1101-SB-2	Stabilizing Bar (Hook)
1	707-L	Adjustable Lateral Support w/11" Arm
3	409	<i>imp</i> ® Rail Clamp

<b>1101-L UPGRADE SYSTEM KIT</b> <b>UPGRADE SYSTEM-CUST WITH 1101-L (LARGE)</b>		
QUANTITY	CATALOG #	DESCRIPTION
1	1101-Case	Storage Case
1	1101-BE	Bariatric Extension
1	1101-D	Anterior Vertical Support (Std)
1	1101-B	Crossover Bar (Std)
1	1101-SB-1	Stabilizing Bar (Tube)
1	1101-SB-2	Stabilizing Bar (Hook)
1	707-L	Adjustable Lateral Support w/11" Arm
3	409	<i>imp</i> ® Rail Clamp



De Mayo Hip Positioner® is a registered trademark of Innovative Medical Products  
 De Mayo Hip Positioner® US Patent # 6,820,621  
 Adjustable Lateral Support™ is a trademark of Innovative Medical Products  
*imp*® Rail Clamp is a trademark of Innovative Medical Products  
 © 2011 Innovative Medical Products  
 All rights reserved



87 Spring Lane, Plainville, CT 06062  
 PH: 860-793-0391 FAX: 860-793-8975  
 Toll Free: PH: 800-467-4944 FAX: 888-229-1452  
 www.impmedical.com  
 Email: info@impmedical.com