

De Mayo Universal Distractor™

New Technique for Femoral Distraction

Improves vision and access to posterior compartment

Reduce procedure time for Unicompartmental, MIS Techniques, Arthroscopies, ACL and TKA surgeries

Requires the De Mayo Knee Positioner® that features a single locking lever for precise control of flexion/extension, tilt & rotation

Patient pad protects the popliteal area

Patent Applied For

Clear field of vision allows easier, unobstructed access



Patent Applied For

Broad distribution of pressure with sterile protector pad that will not rotate or roll

Tension release lever to fine-tune spacing or disengage distraction

Hand control for finite adjustment of distraction

US Patent No. 7,380,299

US Patent No. 7,003,827

A Solution for Tight Spaces

The De Mayo Universal Distractor™ applies finite pressure to the femur outside the knee joint. It allows an unobstructed view into the posterior capsule without the use of lamina spreaders, or a second pair of hands to apply traction. Maximum view allows for easy insertion of instruments, trials, and implants while performing Uni, TKA, ACL and Tibial Plateau surgeries.



De Mayo Universal Distractor™



Surgical Photo courtesy of Robert Kennon, M.D., Waterbury Hospital, Waterbury CT

- ◇ Opens medial and lateral compartments without the use of a “stress” post
- ◇ Surgeon controlled distraction of the compartments reducing potential over stress to the patient’s soft tissue

- ◇ Improves access to the posterior aspect of the joint for removal of osteophytes, resected bone, bone cement, and cautery applications

- ◇ Improves MIS techniques for component insertion and prevents component damage



Surgical Photo courtesy of John Keggi, M.D., Waterbury Hospital, Waterbury CT

Ordering Information

Description	Cat. No.
De Mayo Universal Distractor™ (requires Boot w/Distractor Block)	907
Pressure Protector Pad Set (incl. 294-P Boot Pad, Distractor Pad & Cohesive Bandage)	907-P
De Mayo Knee Positioner® with Universal Distractor™ 30”	803-907-ABD
De Mayo Knee Positioner® with Universal Distractor™ 25”	803-S-907-ABD

Contact IMP for additional ordering options

imp® Innovative
Medical Products, Inc.

A GLOBAL LEADER IN PATIENT POSITIONING

87 Spring Lane, Plainville, CT 06062

PH: 860-793-0391 FAX: 860-793-8975

Toll Free PH: 800-467-4944 FAX: 888-229-1452

www.impmedical.com

De Mayo Universal Distractor™ is a trademark of Innovative Medical Products, Inc.
De Mayo Universal Distractor™ - Patent Applied For
Pressure Protector Pad - Patent Applied For
De Mayo Knee Positioner® is a trademark of Innovative Medical Products, Inc.
© 2010 Innovative Medical Products, Inc.
All rights reserved 02/2010

CAT. No. 907-P
INCLUDES:

STERILE



SINGLE USE



LATEX FREE FOAM