


New

imp Expand-A-Peg™ Locking Peg

the new gold safety standard in pegboard positioning



Uniquely designed to expand in the peg hole!



PATENT APPLIED FOR



Designed for *imp* MORPHBOARD® POSITIONING SYSTEM, while compatible with most other pegboards

- Eliminates any motion of the Peg within the peg hole

- Prevents possibility of loosening intra operatively


- Easy locking mechanism



- Locks in the hole of the pegboard




The most stable peg construct in the industry!




Expand-A-Peg™ Locking Peg

imp's gold-colored, aluminum peg locking system comes in four lengths. A complete positioning system including locking pegs and MorphBoard® can be ordered from *imp*, or components can be purchased individually.



105-1-EP



8", 9", 12", 14" Pegs

Ordering Information

Catalog #	Description
105-EP-8	8" Peg
105-EP-9	9" Peg
105-EP-12	12" Peg
105-EP-14	14" Peg
105-EP	MorphBoard® System Kit with Expand-A-Peg™ Locking Peg
105-1-EP	MorphBoard® Center Board System with Expand-A-Peg™ Locking Peg

imp® **Innovative**
Medical Products, Inc.

87 Spring Lane, PO Box 8028
Plainville, CT 06062

PH: (860) 793-0391 FAX: (866) 459-1805

Toll Free PH: 1-800-467-4944 FAX: 1-888-229-1452

www.impmedical.com