

Sterile **Phase 4 Gel™ Splint**


A UNIQUE STERILE FOREARM TRACTION SPLINT FOR SHOULDER ARTHROSCOPY



- Phase 4 Gel™ Splint distributes pressure and shear forces over a large area. Helps protect against skin breakdown



- Holds patient's arm by adhesion rather than compression as seen in traditional methods



- 3 holes allows for different hook placement to control rotation.

- Pre-formed plastic backbone gives conformity to the forearm/wrist anatomic angle

- Loop and Hook closures hold splint on the patients arm and provides versatility

imp

PATENT APPLIED FOR

Phase 4 Gel™ Splint

- Universal sizing reduces inventory

- Low profile reduces space for storage

- Designed to be used with all traction shoulder positioners



- Can be trimmed for sizing

- Compare to other devices that rely strictly on compression of the arm, wrist/hand or fingers

Ordering Information
Cat. # 509-P4
10 /Case

* ILLUSTRATION ONLY -HOOK(S) NOT INCLUDED



Sterile, Double Peel Pack



Phase 4 Gel™ is a trademark of Innovative Medical Products, Inc.
Phase 4 Gel™ Splint -PATENT APPLIED FOR
Surgical Photos courtesy of Alan Reznik, M.D. New Haven CT
© October 2011 All Rights Reserved

imp **Innovative**
Medical Products, Inc.
87 Spring Lane, PO Box 8028
Plainville, CT 06062
PH: (860) 793-0391 FAX: (866) 459-1805
Toll Free PH:1-800-467-4944 FAX: 1-888-229-1452
www.impmedical.com