



## **Innovative Medical Products Introduces a Better Tightening and More User-Friendly OR Table Accessory Rail Clamp**

*With its handle's special rubber coating and perfect ergonomic shape, the EZ GRIP™ Rail Clamp makes tightening accessory clamps easier than ever*

PLAINVILLE, Conn., (October 19, 2009) – [Innovative Medical Products](#), an industry leader in patient positioning devices and pressure redistribution gel pads, has introduced the [EZ GRIP™](#) -- an operating room table accessory clamp that is easier to tighten, and stay tightened, than other rail clamps on the market today. [The EZ GRIP™](#)'s easy tightening power is made possible by its handle's innovative rubber coating and perfect diameter shape, enabling OR personnel to lock down accessories with more force. In addition, the new clamp comes in a bright green color for quick, easy identification in the OR environment.

The versatile [EZ GRIP™](#) Rail Clamp accepts both round and flat end accessory bars for vertical placement, as well as horizontal extension flat end bars. There is enough space in the "jaws" of the clamp to accomplish these different configurations.

"[The EZ GRIP™](#) offers significant improvements in holding capabilities because of the high torque generated by the extra large comfort grip," said Earl Cole, IMP product development manager. "Other grip handles are either star-shaped or are lever designed, making the mechanism harder to tighten or keep tightened. With [EZ GRIP™](#), accessories will stay locked in place because of its superior gripping power aided by the specially designed and hardened threaded portion of the clamp."

Although the [EZ GRIP™ Rail Clamp](#) is especially designed to accept IMP hip and shoulder positioners, the new non-sterile clamp can also accept most other manufacturers' accessory bars or adapters that have flat or round designs.

### **About Innovative Medical Products**

Established in 1983, [Innovative Medical Products](#) is focused on developing innovative products to benefit and improve efficiency in the operating room and hospital clinics where patient stability and positioning are required. Its products are designed to provide accuracy for the surgeon, save time for operating room personnel, and benefit the overall patient surgical experience. For more information, call: 800-467-4944; or visit: <http://www.impmedical.com>.

###

Media contact: Harry McBrien, Maier-True Communication; 860.677.4581; [harry@maier.com](mailto:harry@maier.com)