


imp Sterilization Case System

for the De Mayo Knee Positioner® and Accessories




De Mayo® Knee Positioner Accessories



De Mayo® Knee Positioner Base Only




L 32" x W 12" x H 6.5" 20 lbs



L 32" x W 12" x H 2.5" 22 lbs



• low profile case with folding pins



• low profile cases comply with hospital weight guidelines

• improves handling efficiency in the OR's sterile field

• conveniently stackable, saves storage space

• helps protect equipment



Holds any combination of *imp* boots and sterile clamps! Components sold separately

Ordering Details

Cat.#	Size & Weight with components
706-KIT	Includes Accessories and Base Plate Case
1210	Base Clamp w/ Fold Down Pins- <i>REQUIRED</i>

Designed to accommodate the De Mayo Knee Positioner® and accessories, including the De Mayo Universal Distractor™, when upgraded to the Base Clamp w/ Fold Down Pins (Cat. No. 1210). Contact *imp* for upgrade information.

imp Innovative
 Medical Products, Inc.
 87 Spring Lane, PO Box 8028
 Plainville, CT 06062
 PH: (860) 793-0391 FAX: (866) 459-1805
 Toll Free PH: 1-800-467-4944 FAX: 1-888-229-1452
 www.impmedical.com