


*Innovative Solutions  
for Today's MIS  
Trauma Challenges*

***imp***<sup>®</sup> *Innovative Medical Products, Inc.*



**IMP...THE OPERATIVE WORD IN PATIENT POSITIONING**

*Innovative Medical Products, Inc.—  
known for its advances in patient  
positioning—offers a line of adjunctive  
orthopedic devices devoted to MIS  
trauma surgery. Working closely with  
key orthopedic surgeons, IMP® has  
brought products to market that not  
only deliver O.R. time and cost savings  
—but promote successful surgical  
outcomes as well.*



## An MIS Approach to Femoral Nailing



The unique, beveled "spade" tip of our 3.0 guide wire eliminates foreign debris in the incision


### Femur Finder®

(Cat. #198)

**Assures accurate guide wire placement for IM fixation and minimizes soft tissue trauma**

- Eliminates the need to "eyeball" the correct femoral entry point
- Simplifies aiming, drilling and directing of guide wire accurately into the femoral canal
- IMP's innovative 3.0 mm guide wire features a specially-designed "spade" tip that eliminates foreign debris in the incision. Use of this proprietary guide wire will facilitate advancement, maximizing performance.
- Does not require direct exposure of the proximal femur
- No need to use an awl to open the femoral canal
- Reduces postoperative morbidity and minimizes complications with hip abductor muscles during rehab

US Patent No. 6,309,396B1




# Tibial Distractor™

(Cat. #602)

**Provides distraction with rotation for axial alignment**

- For use with IM nailing, plating and bone grafting procedures
- Simplifies tibial procedures and eliminates time-consuming fracture table setup
- Autoclavable unit folds down to fit into O.R. sterilizers
- Reduces radiation exposure for surgeons and P.A.s
- Sterile latex-free patient protective foam pads (Cat. #602P)



**Non-invasive boot attachment available**

Patent Applied for



## Knee Positioning Triangles™

(Cat. #393)



**Simple, stable positioning for lower extremity surgery**

- Saves O.R. time...no fracture table setup
- Radiolucent
- Eases muscle tension during fracture reduction
- Latex-free, sterile, disposable pads (Cat. #193-P)
- Autoclavable




## Ducharme FX Alignment Device™

(Cat. #701)

**Facilitates fracture reduction for a closed intramedullary nailing**

- Increases leverage applied during alignment of fragments
- Sliding torque lock is only moving part
- Reduces radiation exposure for surgeons and P.A.s
- Radiolucent
- Autoclavable



# IMP® Can Also Help with Your Knee and Hip Recon Surgeries...

## Knee Positioners to Reduce Surgeon Fatigue and Provide Improved Stabilization

### De Mayo Knee Positioner™ (Cat. #s 803, 803-S)

Single locking lever allows for precise control of flexion/extension, tilt and rotation. Lever also releases boot for testing of ligamentous stability and range of motion. Positioner is easily set up on top of sterile drapes. Available in 30" and 25".




### Universal Knee Positioner™ (Cat. #s 294, 294-S)

Simple adjustment provides superior access to medial and lateral compartments of the knee. Allows for greater adjustability, control and fixation in any desired position. Boot is easily released for testing of ligamentous stability and range of motion. Positioner is set up within sterile field over drapes. Available in 30" and 25".



Distributed By:



## Hip Positioners Ideal for Surgeries Requiring Lateral Decubitus Positioning

### De Mayo Hip Positioner™ (Cat. #1101)

Allows for simple, adjustable positioning. Supports large and obese patients easily. Eliminates pressure on patient's abdomen. Is nurse-friendly, making patient lifting unnecessary during positioning. Cross arm swings away for maximum surgical exposure. Attaches to all O.R. tables.



### Universal Lateral Positioner™ (Cat. #1092)

Provides stable, lateral decubitus positioning. Makes possible a reproducible point of reference for primary and revision arthroplasty. Accessories include: Adjustable Superior Arm (Cat. #7503) Cassette Holder (Cat. #899) Steri-Clamp™ (Cat. #194) Large Model (Cat. #1092L) Dual Compression Foam Pads (Cat. #1092P) Hyperflexion Plate™ (Cat. #298)



**imp**® *Innovative*  
Medical Products, Inc.

87 Spring Lane, Plainville, CT 06062

PH: 860-793-0391 FAX: 860-793-8975

Toll Free: PH: 800-467-4944 FAX: 888-229-1452

[www.impmedical.com](http://www.impmedical.com)

Email: [info@impmedical.com](mailto:info@impmedical.com)