

imp **DEMAYO KNEE
POSITIONER**



Patent applied for

***Single Locking Lever
Allows for Precise
Control of Flexion,
Extension, Tilt and
Rotation***

**IMP...THE OPERATIVE WORD
IN PATIENT POSITIONING**

The DeMayo Knee Positioner by IMP reduces surgeon fatigue and assists in the stabilization of a patient's leg during surgery. The positioner's **single locking lever** allows for precise control of flexion, extension, tilt, and rotation. The same lever releases the boot for testing of ligamentous stability and ROM (Range of Motion) testing. The positioner is easily set up in the sterile field on top of the sterile drapes. Simplicity of use saves valuable time to the O.R. Two sizes are available, 30" and 25" The 25" size will fit into the O.R. Flash Sterilizer. (Cat. #803, 803-S / Short)



Flexion thru extension.



The "flex boot" conforms for a better fit while the low cut design allows for palpation of the maleoli when using an external alignment device. The **Single Lever** controls finite internal / external rotation allowing superior access to the lateral and medial compartments.

Positioner with Single Locking Lever



Controls flexion, extension and boot release.

Single Lever O.R. Table Clamp reduces set up time and securely attaches positioner to O.R. table siderail. Jaws capture sterile drape without perforation. (Cat. #903)



O.R. Table Pads create a well for the positioner. 3 Sectional Pads (Cat. #398)



Sterile Latex-free foam and cohesive wrap secures patient's foot into boot. Packaged 6/Box. (Cat. #294-P)